



Krish Biomedicals

An ISO 9001:2008 Co

Hb/Hct Analyser



Krish Biomedicals

S-35, Pankaj Central Market, I.P. Extension, Patparganj, New Delhi-110092
Tel: +91-11-42648385, Fax: +91-11-42175785, Email: sharadkbpl@yahoo.com

Hb/Hct Analyser

Features

- Compact handy instrument with sturdy design.
- Dual power supply mode with rechargeable batteries
- User friendly functional keys
- 2 line LCD display with convenient view angle
- Automatic self check at start up
- Data storage facility for 200 test results
- Flexibility to use at point of care testing, intensive care units, neonatal wards, external clinics, blood banks
- Maintenance free



*Pre filled ready to use single test
Hb/Hct Analyser HB cuvetts (Reagent cuvetts)*

*Easy to operate with 3 simple steps
1. Open cap of Hb/Hct Analyser Cuvet
2. Add sample, mix and wait for 5 mins.
3. Insert in Hb/ Hct Analyser*

Get result within 5 seconds.

Specification

- * Principle Colorimetry- Cyanmethamoglobin method
- * Light Source: LED
- * Wavelength: 578 nm
- * Path length: 10 mm
- * Sample Volume: 10 μ L
- * Measurement Range: 0-25 g/dL, 3-75 g/dL
- * Analysis Time: Less than 5 second after user request
- * Memory: Upto 200 results
- * Absorbance Range: -0.05 to 2.00
- * Display: Moon light backlit



Krish Biomedicals

S-35, Pankaj Central Market, I.P. Extension, Patparganj, New Delhi-110092
Tel: +91-11-42648385, Fax: +91-11-42175785, Email: sharadkbpl@yahoo.com