



Krish Biomedicals, New Delhi

HEMBRO⁺ Clinical Chemistry Analyser with auto test selection & built in battery Backup

Hembro Plus

The Hembro Plus is an state of art micro controller based analyser designed to perform routine & special Biochemistry test quickly and precisely without the need to select test or setup instrument parameters to run a test. The unique memory chip provided in the reagent kits makes the Hembro Plus ready for the test automatically and stores all the results of the test done.

Instrument Specifications:

Life long LED / Tri-LED Photometer & Silicon photo diode for highly reproducible & reliable results
Results reading in less than 2 Seconds.

Chip Memory to eliminate test selection & programming and to store results of all tests done.

Rs232 port for optional external printer.

Built in Battery Back up sufficient for over 1000 test.

Absorbance Range 0 – 2.5 Abs, Photometric Precision ± 0.002 Abs

Large size Display for easy reading.

Test Available:

HbA1(Ghb)

Hemoglobin

Glucose

Total Bilirubin

Direct Bilirubin

Cholesterol

Triglycerides

Total Protein

Uric Acid

Urea

* Soon Available

- Magnesium
 - Calcium
 - ASO[#]
 - RF[#]
 - CRP[#]
 - HbA_{1c}[#]
 - Lipo Protien(a)[#]
 - Micro Albumin[#]
- [#] Latex Turbidimetry Method



Physical Data:

Dimension : W:228 x D:216 x H:76 (mm)

Printer Paper : Thermal 57mm Wide

Display : 2x16 Chrs Alpha Numeric

Power Supply : Remote 12V 1Amp Adaptor

110 - 220 V AC

Sample Cell: 10 mm Square Glass / Plastic

Up gradation is a continuous process. Specifications and/or appearance are subject to change without any prior notice.